

Hormone Hotline: 228.365.6885



VAGINAL APPLICATOR

VAGINAL APPLICATOR: will help you to deliver the right amount of cream into your vagina. The applicator should only be used in the vagina of adult women.

- 5 GM
- Clear plastic vaginal applicator.
- Designed to be screwed onto the vaginal cream tube.
- Easy and convenient to use.

DIRECTIONS FOR USE:

Remove the cap from the tube of cream/gel & screw the opposite end to the plunger end onto this tube. Squeeze the tube until the applicator barrel is filled with the correct dose ordered. This will push the plunger out. Unscrew/remove the tube of cream/gel. Gently insert the open end of the applicator into the vagina and push the plunger firmly but gently. Keep the plunger firmly pressed down while removing. Wash hands and applicator after use. Use of a sanitary towel is optional to prevent your clothes from getting stained.

TO CLEAN THE REUSABLE APPLICATOR:

- Turn on the hot water to a comfortable temperature.
- Pump the applicator under the running water.
- Take the applicator apart (pull out the plunger from the barrel).
- Wash the applicator parts with soap and water; you may use a Q-tip for removing leftover residue from the tip of the barrel, or better yet-get an eyebrow brush and dedicate it for this very purpose – it works great!
- Rinse well and pat dry.
- Leave disassemble, wrap in a paper towel, and let dry until ready to be reused.

Do not boil the applicator as it will cause it to melt and lose its shape and functionality.

WARNINGS AND PRECAUTIONS:

- Use exactly as instructed by your doctor/pharmacist.
- Seek medical advice if unsure.
- Only for use in adult women.
- If applicator breaks, replace before next use.
- Clean after each use.